

# Lyric® Candidacy

Client name \_\_\_\_\_ DOB: \_\_\_\_\_

Provider \_\_\_\_\_ Date \_\_\_\_\_

## 1. Medical contraindication

(The client is not a candidate for Lyric if the answer is 'yes' to any of the questions below)

Yes No

- Have you ever had radiation therapy to the head or neck?
- Do you have a bleeding disorder (i.e. hemophilia A/B, Von Willebrand, etc.)?
- Do you have chronic ear pain? If yes,  right ear  left ear
- Do you have acute or chronic ear drainage? If yes  right ear  left ear
- Do you have acute or chronic outer, middle, or inner ear disease/disorder? If yes,  right ear  left ear

Is there evidence of:

- Perforated tympanic membrane? If yes,  right ear  left ear
- Malformation or deformity of the ear canal? If yes,  right ear  left ear  
(i.e. interferes with Lyric placement)
- Is there insufficient cognitive ability to understand Lyric device use?

## 2. Additional medical clearance suggested for Lyric candidacy

(Please refer to the appropriate/managing physician if the answer is 'yes' to any of the questions below)

Yes No


- Do you have uncontrolled or insulin-dependent diabetes?
- Do you take any anticoagulant medications (blood-thinners)?
- Do you have a compromised immune system?
- Do you have an implantable medical device (i.e. pacemaker, deep brain stimulation, etc.)?
- Have you had chemotherapy within the last 6 months?
- Is the client under the age of 21?

## 3. Noteworthy considerations

(Extra caution and counseling may be necessary if the answer is 'yes' to any of the questions below)

Yes No

- Do you have a need for regular MRI testing?
- Do you regularly swim under water?
- Do you scuba dive or sky dive?
- Do you have contact dermatitis?
- Is there evidence of exostosis? If yes,  right ear  left ear
- Is there evidence of an osteoma? If yes,  right ear  left ear
- Is there evidence of thin or dry skin in the ear canal? If yes,  right ear  left ear
- Is the ear canal geometry unsuitable (i.e. bulges, vertical step, V-shaped, too narrow)?  
If yes,  right ear  left ear
- Is the audiogram within the Lyric fitting range? If no,  right ear  left ear



Lyric Sizing	Right	Left
Measured insertion depth	mm	mm
Achieved insertion depth	mm	mm
Device size		

**Can the client wear Lyric?**

Yes

Medical clearance required due to: \_\_\_\_\_

Physician name and discipline \_\_\_\_\_

No, due to: \_\_\_\_\_



Lyric, Distributed by Phonak, LLC  
 4520 Weaver Parkway, Warrenton, IL 60555  
 www.lyrichearing.com | 1.800.777.7333 (Option 3) | FAX: 1.866.285.6264

